



Manual HDO - UK/US



# VET-HDO<sup>®</sup>-MONITOR

High Definition Oscillometry

Non Invasive Blood Pressure measurement in dogs and cats



## TECHNICAL DATAS

SAP - DAP - MAP 15 - 300 mmHg  
Heart Rate 550 bpm  
Valve accuracy 15 - 300 mmHg  
Gain - Signal amplification  
Loop function  
PC remote control

## VET-HDO<sup>®</sup>-MONITOR

High Definition Oscillometry

MD PRO



MD / BT PRO



MD Equine



## An important element of:

Routine examinations  
Search for diagnosis  
Control disease / treatment  
Monitoring SVR  
(Stroke Vascular Resistance)  
Monitoring during anesthesia

## SUPPLY ITEMS

VET HDO MD PRO  
VET HDO MD/BT PRO Bluetooth  
for small animals  
VET HDO MD/BT Equine  
✓ BOX  
✓ MDSWIN / HDO Manual  
✓ USB 2.0 Transfer cable

### optional:

BT1000 external Bluetooth-Adapter

✓ Power supply  
INPUT 100 V - 240 V; 0,4 A 47 - 63 Hz;  
Output 5 V 1.75 mA

Battery powered possible  
4 x AA LR6 1,5 V  
Rechargeable Batteries  
4 X AA 1,5 V; min. 2400 mA  
Powerbank 2600 mA

## SUPPLY ITEMS

### Modell:

- ☐ VET HDO MD PRO
- ☐ VET HDO MD PRO / BT
- ☐ VET HDO MD PRO / TAB usb
- ☐ VET HDO MD PRO / TAB BT
- ☐ VET HDO MD Equine / TAB BT

Serial Number: \_ \_ \_ \_ \_

BT1000 Serial Number \_ \_ \_ \_ \_

### Cuffs: Tail measurement is required

- ✓ Standard at VET HDO MD PRO
- C1 - Cat and small dogs ☐ ✓
- D1 - Small dogs ☐ ✓
- D2 - Large dogs ☐ ✓
- H1 - Equine ☐

MDSWIN Software for analysis  
HDO Firmware 2.44  
WIN 8 / WIN 10  
32 / 64 Bit

## WARRANTY

Each System is warranted against manufacturers' defects in workmanship and materials, under requirement use and with proper maintenance for a period of two (2) years after the date of original purchase.

### Date of purchasing HDO

\_ . \_ . \_ . 2 \_ . \_ .

### Distributed by:

\_\_\_\_\_

# TECHNICAL DATA SHEET

<b>Function and operation of each part:</b> Brief explanations of each function
<b>Raw materials:</b> Generic name and their standards of each part mentioned in #1
<b>Electrical Standard / Classification</b>
<b>Manufacturer</b>
<b>Normal rated voltage of mains power supply</b> (optional)
<b>Normal rated voltage of the device</b> HDO MD PRO
<b>Electrical power consumption</b> VET HDO Monitor
<b>Type of Protection</b>
<b>Degree of protection against electrical shock</b> Degrees of protection provided by enclosures (IP Code). Protection of patient against electrical shock classification of the applied part (HUMAN)
<b>Isolation class</b>
Ambient conditions
Storage/transport
<b>System/Block Diagram</b>
Block diagramm
System diagram (Isolations diagram )
<b>Controll</b>
Test and test standard
S + B medVET have entrusted the European Notified Body EUROCAT and METRONIC with identification
<b>In-house test method</b>
Warning for use
<b>Production flow chart</b>
(if there are several manufacturing plants, indicate which plant covers which process)
<b>General</b>
Shut down with Powersupply / Shut down with Battery powered
In the interest of environmental protection, please dispose of the device and batteries in accordance with national or local regulations or environmental protection and recycling



2019



See user manual Hardware HDO- for Software MDSwin- US / UK / DE / IT / ESP / NL  
Rawmaterials\_HDO-US / UK / DE / IT / ESP / NL

S + B medVET GmbH , Neuer Weg 4 64832 Babenhausen, Germany

110-240V AC 1,5 A 50/60 Hz power supply for VET HDO Monitor 6V DC  
(internal batteries: 4x 1,5V AA // or external power supply 6V DC)

5 Watt

Class II

IP 30 Applied part type B without defi-proof Power supply already approved, class 3 internal pump not safety relevant because powered with SELV (safety extra-low voltage)

Isulation class C

EN 60601-1:2008

+10 bis +40 °C 30-75% atm humidity 700-1060 hPa

see: Blockdiagram MD\_UK, Schematic MD\_US

see: Isolationdiagram\_MD, HDO is isolated

CE 0535 o test the electrical safety and the accuracy of data according the relevant standards IEC/EN 60601-1, IEC/EN 1060-1+ 1060-3, Biocompatibility is shown by ISO 10993-1, A risk management has been performed according DIN EN ISO 14971

and their standard related to performance, function and working of device

see: user manual hardware / software

see: Flowchart of production


after 4 minuten in Standby,

Shut down by disconnect the Powersupply

Label according to



## SETUP / DEFAULT





Standby function will be after 3,5 min. active.

The HDO is sleeping but ready and continues to consume electricity!

START


## Activate device

Changes of the basic settings possible via Function menu parameter **FP**



- P0** Adjust the parameters after the first measurement 1 on / 0 off
- P1** Inflates up to 225 mmHg 150 bis 300
- P2** Deflation rate 9 mmHg/sec. 3 bis 18
- P3** Deflates up to 25 mmHg/end of measurement 15 bis 50
- P4** Gain 70 / 100 / 140 / 200 / 280 / 400 / 560 / 640
- P5** Cuff selection C1 / D1 / D2 / H1  
(H1 is only shown by HDO Equine)

## KEYBOARD



START

**Starts the system**

□ There is no data in the memory  
1-6□ There is still data in the memory



**Scroll between the  
Cuffs C1 / D1 / D2**




ENTER

**confirm cuff**

The system is ready to measure.

START

**starts the measurement**



Scroll up

Changes the  
window and enter  
the function menu

Scroll down

Start - Stop - Abort

START

**system wake up, system is booting**


ENTER

one click: confirms the selected cuff  
double click:  
starts the function menu/submenu Parameters

START

**starts measurement**  
In active measurement:  
press START for stop a active measurement  
or finsh / function menu exit

## DISPLAY




- ♥ Flashes - system records
- ⚡ Appears - weak battery
- [L] Automatic cycle (loop) is activated
- [Icon: square with circle] Display of average data
- 1 2 3 There are at least 1/2/3 or more measurements saved in the memory
- [Icon: square with circle] With 3 and more values an average will be calculated **FR**

- FH** Function History = Displays data in the memory
- FR** Function Average = Displays the averages
- FC** Function Clear = Deletes all data
- FE** Function Erase = Deletes single data of your choice
- FL** Function Loop = Set up automatic intervals (Loop)
- FP** Function Parameter = Set up parameters



Up and down between functions

## BASIC SETTINGS




Changes to the basic setting possible via function menu **FP**

## PARAMETER


- P0** Adjust the parameters after the first measurement 1 on / 0 off
- P1** Inflates up to 225 mmHg 150 bis 300
- P2** Deflation rate 9 mmHg/sec. 3 bis 18
- P3** Deflates up to 25 mmHg/end of measurement 15 bis 50
- P4** Gain 70 / 100 / 140 / 200 / 280 / 400 / 560 / 640
- P5** Cuff selection C1 / D1 / D2

## FUNCTION MENU HISTORY



No data in the memory

**1** flashes when the **second** measured value recorded and automatically was saved.



Values and positions are displayed alternatengly

First measured value in the system  
maximum 60 values  
saved

There are 2 values in the memory (history)




Up and down between  
values 1 – 60



back to the function menu




function menu exit



Always the last value measured will be  
displayed  
There are 2 values in the memory


## FUNCTION MENU AVERAGE



Average Value and total number of values  
are displayed alternatengly

Average of 9 measurements

They are average values in the system



There are at least 3 measurements in the  
memory. This is the minimum for an average  
calculation.




Up and down  
between average  
of each position




back to function menu



function menu exit



**FUNCTION MENU CLEAR**




Delete all data

START

starts a new function mode

ENTER

## Flashes




There are 9 measurements in the memory

## Elaborate

## Flashes

-1-

## Flashes



▼

Deletes all data


1

## Flashes

START

```
function menu exit
```

## FUNCTION ERASE



Deletes single data of your choice

ENTER

### Submenu of the parameter FE

Fourth Value

ENTER

## Changes single-reading and position



Up / Down between readings



ENTER


### Confirms single-reading

Display flashes



## Deletes reading

(Now **FA** will display the average of the remaining measurements)



ENTER


back to the display or

START

```
function menu exit
```




## FUNCTION MENU PARAMETER



Change Pre-settings

**ENTER** submenu of parameter FP

PO - automatic AUT  
0 - automatic AUT off  
1 - automatic AUT on



<b>P0</b> Automatic AUT 0 = off 1 = on	<b>P1</b> inflate up to 150-300 mmHg	<b>P2</b> deflate up to 10-50 mmHg
<b>P3</b> deflation rate 3-18 mmHg / sec.	<b>P4</b> Scan with a gain of 70-640 Gain	<b>P5</b> Cuff C1/D1/D2

flashes:  
the value can be changed.



**Up/Down  
between values**




**ENTER** confirms pre-set value

**ENTER** Value is set/saved


Confirm all 4 parameters with

**START**

function menu exit



## FUNCTION MENU LOOP



Automatic measuring interval 1-9 min.


**ENTER** Input start time „t“  
1-60 min.

Display: loop active

Flashes 0 = no loop activated



**Up/down  
between minutes**



Flashes- Value is set to 5 minutes.  
Every 5 minutes a measurement is  
automatically taken.

**ENTER**

confirm set value

n = interval of number of measurements  
consecutively 1-9



**Up / down between  
the values**




**ENTER**


confirm set value

**START**

function menu exit



## FUNCTION RUT



Calculatory measurement


0 - RUT Messung **AUS**  
1 - RUT Messung **AN**

After switching on, the HDO is preset.  
(see basic setting of the system)

What does **RUT** mean?

**RUT** is an automatic calculatory measurement  
which **is not saved!**


An **RUT** measurement matches the HDO setup of the function parameters Pump up Pmax higher and gain (gain / magnification) of the determined values from the first measurement. The basic setting of the parameters depends on the cuff selection.



Which setting has been selected by the HDO is alternately shown in the display with the measurement. If no change is made, the default setting is optimized. An **RUT** measurement is always the first measurements after the start and is only displayed and not saved!

**RUT** stands for one  
Calculatory measurement

## E-ERROR DISPLAY



Error **E1 - E6**

Measurement is not saved!

- E 1** abort by the user
- E 2** too many artefacts
- E 3** Signal amplitude too small, Gain (magnification like a magnifying glass) is automatically increased by one step
- E 4** Deflation rate is too slow
- E 5** ---
- E 6** cuff error:
  - Wrong cuff selected and connected
  - laid out too loosely
  - Cuff defective
  - Valve defective

**INFO:** The HDO shuts off after 5 seconds if no pressure of 20 mmHg was reached or if it was to  
Inflate for more than 8 seconds choose inflation pressure.

## BASICS OF HDO TECHNOLOGY

HDO High Definition Oscillometry has been one for several years patented and commercially available technology for accurate Measurement of systolic, diastolic and mean arterial pressure.

HDO is an evolution of oscillometry and works with high-speed processors and sensors.

It replaces previous technologies (Doppler, conventional oscillometry eg PetMap, Cardell, V20 or Memoprint).

This was the first time that the previous trade fair restrictions were met be lifted by non-invasive blood pressure monitors.

The measurement accuracy of the new technology was added the most diverse species against invasive measuring systems (Gold Standard) confirmed.

### PWA PULSE WAVE ANALYSIS

By using the MDSWIN software analysis, each measurement can be followed in real time on the PC monitor. These pulse waves are further evaluated at the end of a measurement and the overall result then appears as a graphic on the screen.

Each measurement is in the form of values but also as graphics automatically stored in the patient file created for this purpose and can therefore be used later evaluated further.

This technique, due to its 32bit processor capacity a high-frequency direct analysis of the incoming Signals, which in turn leads to a programming of the electronic valves in real time (several times per microsecond), on the other hand allows the assessment of the individual pulse waves similar to the resolution in a CT.

Artifacts and arrhythmias only rarely disturb a measurement or are now assessable. Due to the simultaneous amplitude representation in real time during the measurement can be additional parameters be visually analyzed.


### MORE INFORMATIONS:

[www.hdo-analyse.com](http://www.hdo-analyse.com)

[www.submedvet.de](http://www.submedvet.de)

Download portal nach Registrierung:

<http://www.vethdo.com>



HDO 02-19-005 UK / US